EN - DISPOSABLE POLYPECTOMY SNARES

FR - ANSE À POLYPECTOMIE JETABLE

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DESECHABLES

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Current Name: IFU-Disposable Polypectomy Snares **Creation Date:** 07/02/2022 **Version** 15

This document is only controlled in electronic format by the QM.

EN - INSTRUCTION FOR USE	3
FR - MODE D'EMPLOI	6
NL - GEBRUIKSAANWIJZING	10
ES - INSTRUCCIONES DE USO	13
DE - GEBRAUCHSANWEISUNG.....	17
PT - INSTRUÇÕES DE USO.....	21
IT - ISTRUZIONI PER L'USO	24
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ	28
AR - العربية - تعليمات الاستخدام.....	32
TR - KULLANIM TALIMATLARI	36
SV - BRUKSANVISNING	39

EN - INSTRUCTION FOR USE

Intended Use

Disposable Polypectomy Snares are intended to be used to remove and or cauterize any kind of polyps in the gastrointestinal tract.

Product Description

Disposable Polypectomy Snares are made of a flexible tube with a loop in different sizes and shapes at the distal end and a handle to expose and retract the loop into the catheter. The handle has a HF connector to be attached to a high frequency (HF) generator and an active cord.

Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

Handling and Operation

Insertion into the Endoscope

The Disposable Polypectomy Snare must be passed through the endoscope operating channel with its snare withdrawn into the tube. Please observe compatibility between the Disposable Polypectomy Snare and the operating channel diameters.

HF active cord Connection

Connect the active cord to electricity connector at the handle and to the HF-generator.

HF Generators

Disposable Polypectomy Snares can be connected to all licensed HF monopolar generators by using the active cord.

Always check compatibility between HF generator and active cord. In doubt contact the manufacturer of generator or your local distributor.

The devices are compatible with the following ERBE and OLYMPUS generators:

- ERBE E-Serie and T-Serie

- ERBE ICC
- OLYMPUS HF 120/300

Please refer to the respective manufacturer's manual.

If HF from other manufacturers are used, check their compatibility and follow their instructions. Applicable limit values:

- Max. Power: 300W
- Max. HF voltage: 1600 V_{peak}

Polypectomy snares are composed by a catheter and should never be manipulated with its catheter winded as this may damage the device and make its usage impossible!

Make sure to have on hand a grasping forceps in case the snare get stuck in the tissue and needs to be removed by other means.

Procedure

Clear the area from mucus or other fluids around the polyp as it may decrease the current density reducing the cutting power.

- Place the distal tip of the tube right in front of the polyp to be cut.
- By pushing the handle spool expose the snare fully out of the tube and capture the polyp into the open snare.
- Pull back the spool until the snare closes on the polyp stalk and stays tighten around the polyp pedicle.
- Apply HF-electric current on the snare (Please observe HF generator manufacturer's instructions!).
- Cut the polyp by pulling the snare completely into the tube.
- Stop the HF-electric current and remove carefully the Disposable Polypectomy Snare out of the endoscope.

Do not touch other tissue areas with the snare while using HF-electric current (high hurt risk for the patient)!

Contraindications

- Missing accordance and failing cooperation of the patients
- High – florid inflammation of the colon, toxic megacolon
- Peritonitis, acute abdomen
- Intestine perforation
- Ileus
- Sepsis
- Severe cardio – pulmonary diseases
- Insufficient preparation (cleaning)

Possible complications

Possible complications associated with current application include perforation and bleeding.

Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

Cleaning and Sterilization

This device is designed and warranted for single use!

Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

Serious Incident, Complaints and repairs

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

Declarations

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.
We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

Product disposal

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Legal foundation

The Law of the European Union applied.

Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

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