**EN - DISPOSABLE BIOPSY FORCEPS**

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**NL - WEGWERP-BIOPSIETANG**

**ES - FORCEPS DE BIOPSIA DESECHABLE**

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**Current Name:** IFU-DisposableBiopsyForceps

**Creation Date:** 2024-06-25 **Version** 19

*This document is only controlled in electronic format by the QM.*

EN - Catalog number		
FR - Numéro de catalogue		
NL - Catalogusnummer	U120C	U230C
ES - Número de catálogo	U160	U230SC
DE - Katalog-Nummer	U160C	U422-180
PT - Número do catalogo	U160PC	U422-180S
IT - Numero di catalogo	U160PCS	U422-230
EL - Αριθμός καταλόγου	U160SC	U422-230LC
AR - العربية- رقم الكتالوج	U230	U422-230S
TR - Katalog numarası		
SV - Katalognummer		

EN - INSTRUCTION FOR USE .....	3
FR - MODE D'EMPLOI .....	6
NL - GEBRUIKSAANWIJZING .....	9
ES - INSTRUCCIONES DE USO .....	13
DE - GEBRAUCHSANWEISUNG.....	16
PT - INSTRUÇÕES DE USO.....	19
IT - ISTRUZIONI PER L'USO .....	23
EL - ΟΔΗΓΊΕΣ ΧΡΉΣΗΣ .....	26
AR العربية- تعليمات الاستخدام.....	30
TR - KULLANIM TALIMATLARI .....	33
SV - BRUKSANVISNING .....	36

# EN - INSTRUCTION FOR USE

## 1. Intended Use

Disposable Biopsy Forceps are used for taking out biopsy specimens (tissues samples) from the human body, through endoscopes.

## 2. Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

## 3. Contraindication

All the usual contraindications for the use of disposable biopsy forceps apply to the use of G-Flex disposable biopsy forceps. These include:

### **Related to the cardiovascular system:**

Shock, Hypoxia, Comorbidity (e.g. severe cardiopulmonary diseases and decompensation), Coagulopathy, Uncontrollable hemorrhagic diathesis, Acute myocardial infarction, Tissue samples from vessel injury in ulceration, Use on the central cardiovascular system.

### **Related to the gastrointestinal system:**

Peritonitis, Fragility of the intestinal wall (e.g. ulcerative colitis, diverticulitis, colitis ulcerosa, toxic megacolon), Fulminant colitis, Ileus, Not prepped bowel, Recent gastrointestinal anastomosis, Acute abdomen.

### **General contraindications:**

Non fasted patient, Severe respiratory insufficiency, Acute perforation of the tracts; Sepsis, Pregnancy.

Moreover, an analysis of the risk/benefit ratio should be established by endoscopist and cardiologist for patients using anticoagulants or anti-aggregants.

## 4. Possible complication

Possible complications are perforation, bleeding, difficulties to remove/insert, mechanical issues, bent device issues, breakage issues.

## 5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## 6. Product Description

The device is composed by a jaw (in different shapes and sizes), a spiral spring hose (in different lengths and diameters) and an operating handle.

## 7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

Unwind the disposable endoscopic forceps and/or retrievers, manipulate backward - forward smoothly the plastic spool of the handle to verify the opening / closing of the cups.

If the package present any damaged, the sterility of the device is not guarantee.

## 8. Handling and Operation

The biopsy forceps must be passed through the endoscope operating channel with its cups closed. Please observe compatibility between the biopsy forceps and the endoscope operating channel diameter.

Unwind the disposable biopsy forceps, manipulate backward - forward the plastic spool of the handle to verify the opening / closing of the cups.

## 9. Procedure

- Approach the biopsy forceps near to the tissue (lesion) to be taken from the patient's body.
- Open the cups by pushing the plastic spool at the handle.
- Place the forceps on the tissue and take the sample gently by pulling back the plastic spool (this manipulation will close the forceps cups).
- Keep the forceps cups closed and remove biopsy forceps out of the endoscope carefully.
- This tissue specimen taken will be ready for further histological examination.

## 10. Precautions/ Warnings

Disposable Biopsy Forceps composed by a spiral and should never be manipulated with their spiral winded as this may damage the device and make its usage impossible.

## 11. Cleaning and Sterilization

This device is designed and warranted for single use!

## 12. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!  
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

### **13. Risk in case of reuse**

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

### **14. Serious Incident, Complaints and repairs**

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

### **15. Declarations**

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

### **16. Product disposal**

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

### **17. Legal foundation**

The Law of the European Union applied.

## 18. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday; 9am to 5pm (MET)

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