**EN - EXTRACTION BASKETS & LITHOTRIPSY SYSTEM**

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


**Current Name:** IFU-Basket & Lithotripsy System

**Creation Date:** 07/02/2022


**Version:** 4

*This document is only controlled in electronic format by the QM.*

EN - INSTRUCTION FOR USE .....	6
FR - MODE D'EMPLOI .....	13
NL - GEBRUIKSAANWIJZING .....	20
ES - INSTRUCCIONES DE USO .....	27
DE - GEBRAUCHSANWEISUNG.....	34
PT - INSTRUÇÕES DE USO.....	41
IT - ISTRUZIONI PER L'USO .....	48
EL - ΟΔΗΓΪΕΣ ΧΡΪΗΣΗΣ .....	55
AR العربية - تعليمات الاستخدام.....	63
TR - KULLANIM TALIMATLARI .....	70
SV - BRUKSANVISNING .....	76

# EN - INSTRUCTION FOR USE

## 1. Intended Use

Extraction Basket & Lithotripsy System are intended to catch and remove stones or to crush stones in the biliary and pancreatic ducts of the human body.

## 2. Intended Population

Our device is intended for both adult and paediatric patient according to the physician recommendation taking into account the contraindications.

## 3. Contraindications

There are no contraindications currently for this device family.

## 4. Possible Complications

Possible complications associated with perforation, bleeding, pancreatitis, cholangitis and cholecystitis.

## 5. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

## 6. Product Description

The **extraction basket** is made of a basket at the distal end and a handle at the proximal end. Some models have a removable handle (see page 3)

1. Basket
2. Catheter
3. Luer-lock connection
4. Locking system of the removable handle
5. Removable handle

The **spiral** is made of:

6. Rigid shaft
7. Spiral

The **lithotripsy handle** is made of a winding cylinder with a locking system fixed to a wheel, an aluminum body with an opening for the spiral and a handle (see page

- 4). The lithotripsy handle is delivered non-sterile and is made of the following parts:
8. Spiral entrance
9. Winding cylinder with orifice for tracking wire
10. Locking system
11. Wheel
12. Handle

## 7. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential. Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package presents any damage, the sterility of the device is not guaranteed.

## 8. Handling and Operation

The device has to be withdrawn completely into the Teflon tube to be introduced into the endoscope operating channel. Please observe compatibility between the device and the endoscope operating channel diameters.

## 9. Procedure

### Extraction of Stone

- Once the stone is located by usual means proceed as following:
- Place the Extraction Basket's distal tip above the stone to be removed.
- Push the basket completely out of the catheter and pull the basket down to the stone.
- If necessary, rotate the basket or move the catheter up and down in short and quick moves in order to facilitate to place the basket wires around the stone.
- Pull the handle back until the stone is fixed in the basket and pull the device toward the endoscope until the stone is extracted.

Excessive traction should be avoided to prevent damage the endoscope.

### Emergency lithotripsy

- In case a stone gets stuck inside of an Extraction Basket proceed as following:
- With the help of a cutting plier, cut the catheter and tracking wire of the basket as close as possible to the handle of the device, just before the shrinking tube.
- Dispose the handle and pull gently the catheter over the cut wire 1-2cm to ensure the wire is protected and can not damage the endoscope.
- While keeping both tracking wire and catheter inside of the patient, gently remove the endoscope. Ensure not to move the catheter up or down while removing the endoscope.
- Once the endoscope is completely removed, remove the plastic catheter and replace it by the short and thick metallic spiral (Ref. GF641).
- Once the spiral reaches the stone, push the tracking wire through the lithotripsy handle spiral entrance until a slight tension can be felt.
- Lock the lithotripsy handle at the locking system.
- Gently begin to rotate the wheel until the stone gets crushed. When applying rotation on the wheel, allow some delays so that the mechanical force is perfectly transferred to the distal end of the spiral before repeating a new rotation ( at least 40 seconds)
- Remove the device from the patient.

It is essential to keep the spiral as straight as possible to avoid friction between the tracking wire and the spiral. Any friction could lead to break the tracking wire and cause the procedure to fail.

Do not rotate the wheel quickly as this may lead to break the tracking wire and cause the procedure to fail.

### **Lithotripsy through the working channel (for baskets GF264GW, GF265GW, GF266GW, GF267GW only)**

#### **How to use the lithotripsy basket with a guide wire**

This model of basket can be introduced into the bile duct directly or over a previously positioned guide wire. To do so, please follow the instructions below:

- Loosen the sealing cap at the Y-adapter.
- Remove the basket pulling it from its handle (without unscrewing it) from the distal outer teflon tube.
- Insert the teflon tube over the guide wire through the working channel of the endoscope until the metallic tip of the tube reaches a position behind the stone to be crushed.
- With the teflon tube in position, remove the guide wire from the proximal end of the teflon tube.
- With two fingers, squeeze the basket and introduce it through the proximal end of the teflon tube (see page 5) until the rigid part of the wire along with the handle reach the teflon tube.
- Slowly expose the basket and catch the aimed stone to be crushed.

#### **Changing the Teflon tube for the metal spiral**

For crushing the stone it is necessary to change the Teflon tube for the fitted metal spiral. These action makes a stable counter-pressure during the crushing procedure possible.

- Loosen the sealing cap at the Y-adapter.
- Remove the Teflon tube completely and carefully out of the endoscope, sliding it over the traction cable.
- Insert the metal spiral's distal end by sliding it over the traction cable into the operating channel until it reaches the basket with the stone inside.

#### **Handle Assembling (See page 4)**

- Keep the locking system in the open position
- Push the traction cable tip through the spiral entrance to the winding cylinder orifice.
- Push the handle over the traction cable until it couples to the spiral tip and insert the spiral tip in the spiral entrance until resistance is felt (excessive force should be avoided).
- Push the spiral straighten to the basket.
- To release the tracking wire, place the locking system in the close position.

## Crushing the Stone

- Turn the manual wheel at the handle, until the stone is crushed. When applying rotation on the wheel, allow some delays so that the mechanical force is perfectly transferred to the distal end of the spiral before repeating a new rotation ( at least 40 seconds)

A lithotripsy has no guarantee of total success! If the stone to be crushed is very hard the basket might be torn because of the force applied on it. In this case the stone along with the torn basket should be surgically removed.

The lithotripsy should be interrupted immediately if the stone crushing could not be achieved.

## 10. Precautions/Warnings

The device has to be withdrawn completely into the Teflon tube to be introduced into the endoscope operating channel.

Please observe compatibility between the device and the endoscope operating channel diameters

## 11. Cleaning and Sterilization

Before the first use and prior to all subsequent uses, the device must be cleaned, disinfected and sterilised. Effective cleaning and disinfection are essential for an effective sterilisation.

The reprocessing of reusable endoscopic equipment must remain under the responsibility of individuals with thorough training in reprocessing methods and be aware of the following points:

- Local hospital's reprocessing guidelines
- National health and safety guidelines
- The instructions in this instruction for use
- Additional national specifications related to effective inactivation of prions (i.e: use of highly alkaline cleaning agent and prolonged sterilisation time)
- That the methods used to clean, disinfected and sterilise are adequately validated for this type of instruments
- That the equipment used is regularly maintained and checked
- That the validated parameters are maintained during each reprocessing cycles

The reprocessing of the instrument must be initiated immediately after use, since dried organic residue can impact effective cleaning, disinfection and sterilisation.

### Reprocessing methods

A mechanical cleaning/disinfection process should be used whenever possible over manual cleaning.

Therefore, one of the following 2 options must be used for the reprocessing:

1. Perform the pre-treatment and manual cleaning before the mechanical cleaning/disinfection followed by the sterilization (3 steps).
2. In case a mechanical cleaning/disinfection is not available the device must then be pre-treated and manually cleaned followed by a sterilization (2 steps).

## Material Resistance

Ensure that the cleaning and disinfection agent do not contain the following:

- Aldehydes (this will fix blood-based soiling)
- Strong organic, mineral or oxidising acids (lowest permissible pH value 5.5)
- Strong alkalis (highest permissible pH 11, neutral / enzymatic or slightly alkaline cleaners are recommended)
- Organic solvents (e.g. ether, ketones, benzines), fluorinated alcohols
- Oxidising agents (e.g. hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic / halogenated hydrocarbons
- Formamide
- Trichloroethylene / perchloroethylene

## Pre-treatment / Manual Cleaning

1. Immersed the instrument in cleaning solution within maximum 2 hours after use
2. Remove any soiling visible on the jaws using a clean soft brush (or a clean, soft and lint-free cloth) which is used solely for this purpose. Never use metal brushes or wire wool.
3. Open and close the jaws, immersed in cleaning solution at least 3 times
4. Immediately place the instrument wined with a diameter not less than 20 cm into an ultrasound bath filled with cleaning solution for 30 minutes (time and concentration may vary according to the solution manufacturer's specifications. Cleaning method has been validated with ultrasound bath using 0.5% neodisher® Mediclean forte (Dr. Weigert GmbH) at 35 kHz)
5. Rinse the instrument for at least 1 minute under running water (temperature < 35°C/85°F)

## Mechanical Cleaning / Disinfection

Proceed as followed:

1. Place the instrument wined with a diameter not less than 20 cm into the cleaning and disinfection equipment that fulfils the following requirement:
  - Equipment is effective and suitable for the cleaning and disinfection of medical instruments (i.e: FDA approved, CE marked)
  - Only sterile or low-germ (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units /ml) water is used (e.g. purified water/highly purified water)
  - The air used for drying is filtered (oil free, low-germ and low-particle)
  - The device is maintained and checked regularly
  - The cleaning and disinfection solutions used are compatible with the instrument (See "Material Resistance")
  - The concentration, temperatures and exposure times specified by the cleaning / disinfectant agent manufacturer is followed at all time, this includes the rinsing specifications
2. Start the cleaning and disinfection program:
  - The programme must be suitable for the instrument

- Thermal disinfection programme should be used (90°C / 194°F for at least 5 minutes)
  - The program contains at least 3 rinsing cycles (depleting steps) after cleaning (or neutralization if applied) or conductance control recommended to effectively prevent detergent residues
3. Remove the instruments from the equipment when the program is finished
  4. Check the functionality of the device by opening and closing the jaws several times
  5. Check the device for signs of corrosion, damaged surfaces, cracks, soiling or discolouration
  6. Remove any damaged device from use
  7. Any soiled device must be cleaned and disinfected again
  8. Package the device into a single-use sterilisation pack or a sterilisation container (suitable for steam sterilisation of at least 134°C / 270 °F) as soon as possible after removal from the equipment

## **Sterilisation**

Only steam sterilisation must be used.

Flash sterilisation, hot air sterilisation, radiation sterilisation, formaldehyde, ethylene oxide sterilisation or plasma sterilisation are not permitted.

The following aspects shall be taken into account:

- Fractionated vacuum method (at least 3 vacuum steps and with adequate product drying not less than 20 minutes)
- Steam steriliser compliant with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ and product-specific performance assessment (PQ))
- Maximum sterilisation temperature 134°C / 273°F plus tolerance as defined in DIN EN ISO 17665
- Sterilisation time (exposure time at sterilisation temperature) at least 4 min at 134°C/270°F).
- For prion inactivation (depending on the local recommendation) at least 18 minutes at 134°C/273°F

## **12. Storage**

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!  
Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

### **13. Risk in case of reuse**

The Disposable extraction basket is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

### **14. Serious Incident, Complaints and repairs**

In case of a serious incident, please report to the manufacturer and to the competent authority of the European Member State where the user is established. Shall the incident occur outside the Europe, please report to the manufacturer and follow the local applicable regulation.

In health protection of our employees only disinfected or sterilised instruments will be accepted for analysis or repairing. In any case disinfection or sterilisation dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

### **15. Declarations**

We, hereby, declare that the medical device neither contain nor is produced with tissue of animal origin.

We, hereby, declare that the medical device neither contain nor is produced with tissue or extracted material of human origin.

We, hereby, declare that the medical device do not contain medicinal substance.

We, hereby, declare that the medical device do not contain human blood derivative.

We, hereby, declare that the medical device do not contain phthalates or other products that could leak from the device that may be carcinogenic, mutagenic or toxic.

We, hereby, declare that the medical device is not composed of substances or combination of substances intended to be introduced into the human body.

### **16. Product disposal**

After use, these products may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

### **17. Legal foundation**

The Law of the European Union applied.

### **18. Product support**

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

**Times available:** Monday to Friday; 9am to 5pm (MET)

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